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PCU vario 40

Translation of the original instructions for use

1. Device description

The PCU vario is a light polymerization device for the curing of components which have been created by means of 3D printing. The following device properties are used in the curing process:

Intended use:
The PCU vario may only be used for curing of resins and lacquers which are crosslinking in a wavelength range between 375 and 450 nm.
The following device properties are used in the curing process:

- LEDs must not be covered.
- Use only in dry environments.
- Operate only with closed polymerisation chamber.
- The needle valve may only be adjusted by trained maintenance personnel.
- The unit may only be used with nitrogen (N_2) as the inert gas.
- A suitable coupling socket must be used on the N_2 connection.
- The current draw of the device connected to the USB port must not exceed 500 mA.
- The maximum permissible length of the connection cable is 3 m.
- Do not lift the unit by the handle of the lid.

Device properties:
In the curing process, the following device properties are applied:

a) LED technology
The advantages of the LED technology are that a permanent light of one wavelength is emitted. The permanently emitted light radiation can be captured by a standard UV sensor and thus allows for the measurement and the observation of the light output in the exposure chamber.

b) Electronic control (micro-controller)
The use of a micro-controller allows for the control and reading of the current, pressure and temperature sensors.

c) Simple control concept
The device is operated by the removable rotary-push-button. The corresponding option is selected by turning the button. The selected option is operated by pressing. Longer pressing opens the parent and sub-folder menu. Both elements can be cleaned with the appropriate cleaning agents after removal.

d) Curing environment
The components to be cured can be cured under different environments. The environment to be set depends on the application and can be found in table 1.

Table 1: Curing environment for different uses

Environment	Use	PCU vario 40
Nitrogen	Curing of medical devices	✓
Ambient air	Curing of laboratory products	✓

2. Technical data

Weight of the overall package	26 kg
Weight of the device	23 kg
Dimensions of the overall package	H 535 x W 740 x D 540 mm
Dimensions of the post-exposure module	H 165 x W 330 x D 540 mm
Dimensions of the post-exposure module with control module	H 210 x W 470 x D 540 mm
Polymerization chamber	H 100 x W 270 x D 200 mm
Usable volume of the polymerisation chamber	H 60 x B 250 x T 195 mm
Voltage	110–240 V AC
Nominal frequency	50 / 60 Hz
Max. Power consumption	1200 W
LED wavelength	375 nm, 405 nm, 450 nm
Light intensity	max. 300 mW/cm² per side
Ambient temperature	Nitrogen, ambient air
Max. altitude	2000 m
Ambient temperature	10–30 °C max. rel. Luftfeuchte 80 %

Fluctuations in the mains supply voltage up to ± 10 % of the nominal voltage transient overvoltages up to the values of overvoltage category II temporary overvoltages.

Permissible degree of soiling	2
Fuse	16 A slow-blow fuse specified to UL 284-14

3. Safety instructions


Attention! Please read these notes carefully before the connection and commissioning of the device. The operational safety and functioning of the device can only be ensured if both the general safety and accident prevention regulations as specified by the legislator and the safety instructions in the operating instructions are observed.

Disconnect the mains plug before cleaning and maintenance work. Risk of electric shock.

In case of long and intensive use or leakage of the inert gas, the oxygen might be displaced. The unit may therefore only be used in well-ventilated environments.

Do not look into the lamp for an extended period of time while it is operating. It can be harmful to the eyes.

The unit heats up during operation. Burn hazard!

• Warning notices and other labels on the appliance must not be removed.

1. The device must only be used in accordance with the present operating instructions. We do not assume any liability for damages which result from improper use and/or faulty operation.

2. Only operate the device at a socket with a protective conductor. Never touch the plug with wet hands.

3. Use only the supplied power cable or a power cable that has the same technical specifications and the necessary certifications for use in your country.

4. The device must be placed on a stable and level surface which is sufficiently stable and loadable for the operating.

5. This device is equipped with LEDs of Risk Group 2 according to IEC 62471 standards. This classification exceeds the exempt group, indicating a potential hazard to eyes and skin. The risk to the viewer depends on how the device is installed and used.

To minimize risks, the device is equipped with a drawer switch that automatically turns off the LEDs when the drawer is opened during proper operation. However, it is still recommended not to look directly at the LEDs and to maintain a safe distance.

6. Do not insert any objects into the device apart from the polymerization chamber.

7. Do not lock the operating and safety elements of the device.

8. Please check the device to ensure its proper condition and operating safety before each operation. If its condition is not flawless, the device must not be used and must be labelled accordingly.

9. Before the cleaning and maintenance of the device or before replacing parts, the power plug must be pulled out.

10. Signs and stickers must always be kept in an easy-to-read condition and must not be removed.

11. The opening of the device and repairs may only be performed by authorized professionals.

12. Only use nitrogen as protective gas. Do not use flammable or toxic gases under any circumstances. Please observe the general safety and accident prevention regulations when handling nitrogen.

13. Only use approved gas pressure bottles.

14. When connecting the gas pressure bottle to the device, be sure to check the tightness of the feed tube. Please make sure that the gas pressure bottle cannot tip over.

15. Never connect the gas bottle without using a pressure reducer; the device may be operated with a maximum pressure of 5 bar.

16. Please close the gas supply, when the device is not used. Do not operate the device unattended.

17. Be sure to provide for sufficient ventilation in the work area.

18. Do not keep any highly inflammable substances within the immediate vicinity of the device.

19. Authorized operators: The operator of the device must make the operating instructions accessible to the operator and must also make sure that s/he has read and understood them. Only then may the operator put the device into operation.

20. To prevent the ingress of water into the device (e. g. splash water), the device should be installed in a dry environment.

21. Please disconnect the device from the mains when it is not used.

22. Do not operate the device if the housing of the device is not closed and screwed on.

23. Unauthorized modifications and changes are inadmissible for safety reasons.

24. The prescribed operating and maintenance conditions specified in these operating instructions must be complied with. Please observe the general accident prevention regulations when working with the PCU vario.

25. Wear gloves when handling uncured components or plastics.

4. Commissioning

Please ensure before the commissioning that no damage has occurred during transport. Check the housing for external damages and control the glass plate in the polymerization chamber. The drawer is secured for transport reasons. The transport locks must be removed before the commissioning of the device. Attach the control module to the designated mounting bracket using the screws supplied. Then connect the operating module and the post-curing module using the included bus cable. You can use up to three post-curing modules with one control module. Please refer to the section on using multiple post-curing modules. Please ensure that the indicated voltage on the type plate corresponds to the voltage of the power source. Connect the power cord to the appliance inlet and connect the unit to the mains using the power cable provided. Attach the protective gas supply to the connection of the pressure reducer of the gas pressure bottle with the use of the hose clamp. Establish a connection to the device with the quick-lock coupling.

Scope of delivery:

1x PCU vario 40
1x Accessory package consisting of:
3x Power cables
Only connect the cable with the plug that matches the local socket system.
1x Operating instructions
1x Gas connection line kit including quick coupling and 2 hose clamps
1x Modbus cable
1x Holder for control module with 4 screws

5. Information on the use of the protective gas


The use of the protective gas allows for the completion of the polymerization in an atmosphere without oxygen. The protective gas replaces the air in the polymerization chamber and thus reduces the oxygen content in the polymerization space of the PCU vario, which leads to the hardening of components without an inhibition layer. The protective gas supply is automatically controlled when the polymerization process starts. The flow rate is preset; the adjustment can be changed at the control needle valve at the back of the device (please observe the safety instructions).
The outlet pressure of the bottle can be mounted to the pressure reducer of the gas pressure bottle so that the flow is changed without changing the control needle valve. The gas flow is factory-set The use of a flow meter is recommended for monitoring the gas flow during operation. (Please refer to the instructions for use of the relevant material for exact information on type and flow rate.) The connection pressure must be set to 5 bars.

6. Functional elements

Overview

Rear view

7.1. Operation: Starting the device

The device is connected to the socket using the mains cable supplied and switched on using the main switch on the back.

After switching on the appliance, the user interface opens. The programme sequence can be set using the rotary/push button.

7.2. Functions of the rotary-push-button

Selection


Confirmation

7.3. Operation

Current home screen


After turning on the device, the home screen is displayed for a few seconds.


The display then automatically switches to the "Curing Profiles" menu.

Navigation	
Curing Profiles	**Settings**
Profil 1	Import Curing Profiles
Profil 2	Get New Profiles
Profil 3	Set Time/Date
	Update Curing Device
	Update Vacuum pump
	Update UI Bootloader
	Info
	Display Dimming


The menu structure of the
PCU vario 40 is organized as followst:

PCU vario 40

Translation of the original instructions for use

The current time is displayed on the left side of the header. There are two icons on the right.

This icon takes you directly to the "Navigation" menu.


This icon opens the parent menu.

The footer displays which post-curing module the control unit is currently communicating with.
A single control unit can control up to three post-curing modules.

Curing Profiles

The following section provides a more detailed explanation of the individual menu items.

Turn the control knob to select a program from the list. Press the control knob to confirm your selection.


If the„ ▶" icon is highlighted, press the control button to start the selected exposure program.

If the "X" icon is selected, you can cancel the currently running program by pressing the control button.


After the program has completed successfully, the following message appears

Curing finished
Press key to continue

on a green background

Settings

This menu allows you to access the submenus for settings.

The menu contains 8 submenus which are described individually below.

- Import Curing Profiles
- Get new profiles
- Set Time/Date
- Update Curing Device
- Update Vacuum pump
- Update UI Bootloader
- Info
- Display Dimming

Import Curing Profiles


Importing Material Profiles from a USB Drive.

Insert a USB drive containing the file "programs_vario.json" into the USB port on the back of the control unit. Then press the control button to start the import process.

Do not remove the USB flash drive until the import is complete. Removing it prematurely may result in an incomplete or corrupted data transfer.

Note: When importing a new program list, the existing program list will be completely overwritten. It can only be restored by re-importing the corresponding JSON file.


Get new profiles

After scanning the QR code shown with a cell phone, provided there is an Internet connection, a web page will open from which curing profiles for the PCU vario 40 can be downloaded.

Set Time/Date


Setting the Date and Time

Turn the control knob until the desired setting is highlighted, then press the control knob. You can now adjust the selected value. Press the control knob again to save the setting.

Update Curing Device


Updating the Post-Curing Unit's Firmware

If necessary, Dreve provides a firmware file and instructions for performing the update. Firmware updates can introduce new features.

Update Vacuum Pump


This feature is not available for the PCU vario 40. It is required exclusively for the PCU vario 160.

Update UI Bootloader


Updating the Bootloader


If necessary, Dreve will provide the required file and instructions for performing the update.


Info

This menu displays information about the software version and the device's serial number.


Display dimming

Adjusting the Display Brightness

7.4. Placing the components

The components are placed on the glass plate in the drawer of the polymerisation chamber.

7.5. Use of Multiple Post-Curing Modules

A single control module can control up to three PCU vario 40 post-curing units. The control module is connected via a bus cable to the post-curing unit on which it is mounted.

Up to two additional post-curing units can be stacked on top of the first post-curing unit. The post-curing units are connected to one another via bus cables. The order of connection is arbitrary.

To ensure correct addressing of the post-curing units in the bus system, a unique device address must be set on the coding switch of each post-curing unit. If two or more post-curing units are configured with the same address, the bus system will not be able to recognize them properly.

The following figure shows the correct setting of the coding switches for assigning device addresses.

Device 1

Device 2

Device 3

8. Maintenance

Maintenance work to be performed by the user

Maintenance unit	Maintenance frequency
Clean the glass panes in the bottom with isopropyl alcohol	daily*
Before commissioning check whether all connections of the device are positioned correctly: • gas connection • mains connection	daily*
Cleaning of the device with a mild detergent	weekly*
Check safety functions at regular intervals:	weekly*

Maintenance work to be performed by the service commissioner

Maintenance unit	Maintenance frequency
Check whether the gas valve works reasonably	weekly*
Check the lid closing function	weekly*

* Recommendation, actual need depends on use intensity and other requirements.

9. Error Codes


The various process parameters are measured by the sensors installed in the PCU vario.

10. Troubleshooting plan

Code	Defect	Cause	Troubleshooting
1	Error Overtemperature	Overtemperature	Let the device cool down
2	Error LED Driver Communication	LED Driver Communication Error	The device must be inspected by a qualified technician
4	Error LED Driver Current/Voltage	Faults in the electronics of the LED drivers	The device must be inspected by a qualified technician
8	Error Pressure Sensor Communication	Error in data communication with the pressure sensor	The device must be inspected by a qualified technician
16	Error Pressure Sensor CRC	Pressure sensor is providing incorrect data	The device must be inspected by a qualified technician
32	Error Modbus Communication	Bus cable not connected properly	Checking the Connector Connections on the Bus Cables
		Coding switch not set correctly	Checking the settings of the coding switches on each device

64	Error Overcurrent	Excessive current draw detected	The device must be inspected by a qualified technician
128	Error Nitrogen Flow	The nitrogen valve on the gas cylinder is not open. Caution: The device may be operated at a maximum pressure of 5 bar (70 psi)	Open the nitrogen valve
	The program won't start	The drawer isn't closed	Close the drawer
		Drawer limit switch defective	Replace drawer limit switch
	Device is not working	Incorrect mains connection	Check mains connection
		Main switch off	Turn on the main switch
		Fuse is defective	Replace the fuse according to the specifications in the technical data

11. Exploded view

Pos.	Designation & Configuration	Item Code
1	PCU Vario Control Module Housing	D521265
2	PCU Vario Control Module Back Panel	D521266
3	Control Module Circuit Board	52301
4	Rotary Knob	521029
5	PCU LED Rotary Knob Adapter	521044
6	PCU Vario Control Module Display Clip	521593
7	Grayhill Encoder 62AG	521566
8	Display	52300
13	Control Module Sheet Metal Bracket	521191

Pos.	Designation & Configuration	Item Code.
1	LED Circuit Board vario 40	521635
3	Mirrored Sheet Metal Chamber vario 40	521664
4	PCU Heat Sink vario 40	521634
6	PCU Chamber Wall vario 40	521627
7	LED Mirrored Sheet Metal vario 40	521632
8	Mini Angled Screw-In Fitting with O-Ring Seal	5533
9	Snap Switch, Limit Switch	521572
10	Ring Magnet 2.0 x 7, Hole 4	521051
11	VA1245 Telescopic Slide, 200 mm	521684
12	PCU vario Drawer Front Side Panel, Left	521655
13	PCU vario Drawer Front Side Panel, Right	521656
14	PCU vario Drawer Handle, 40	D521654
15	Glass top, PCU vario 40	521630
16	Drawer insert, PCU vario 40	521629
17	Drawer front center panel, PCU vario 40	D521621
18	Drawer box, PCU vario 40	521628
19	Drawer seal, vario 40	521668
20	Drawer backing plate, vario 40	521636
21	Ring Magnet Mount for Drawer	521691
23	Ring Magnet Mount for Drawer	521624
24	Connecting Plate for Vario 40 Driver-Control Board	521663
25	Component Carrier for Vario 40	521662
26	Power Distribution Board for Vario 40	521657
27	Control Board for PCU Vario Post-Curing Module	53200
28	Driver Board for Vario 40	521649
29	MSP-1000 Power Supply	521231
30	Mains Input Jack + Filter Schaffner FN282-6-06	521568
31	Feindrosselventil	50937
32	Flow-Restrictor Valve	521640
34	Solenoid Valve A331-1C2 3/2-Way Solenoid Valve 1/8 NC	5534
35	Straight Plug Connector. (5534) 6-1/8-MSvbl	521695
36	Pressure Sensor MPX5100DP	521075
37	1/8" Elbow with Gasket s2020	50095
38	Plug nipple with G1/8 male thread, 5 mm	521637
39	Cover for Vario 40 post-curing chamber	521625
40	Vario 40 PCU side bracket, right	521626
41	Vario 40 PCU side bracket, left	D521665
42	Wiring harness for driver board-IDC adapter, bottom, Vario 40	D521666
43	PCU housing, Vario 40	D521620
44	PCU base plate, Vario 40	521623
45	Spacer strips for post-curing chamber, Vario 40	521686
46	Unit foot	521230
47	Sanyo Denki 9G DC axial fan, ball bearings, 24 V	521680
48	Fan 60x60x35, 48 V, Sanyo Denki	521681
49	RGB indicator light	508356

12. General information

Repair and maintenance work at the PCU Vario may only be carried out by qualified personnel. The power supply to the device must disconnected in any case.
The company Dreve will be responsible for the safety, reliability and performance of the device, if

- the device is only used for the appropriate application purpose.
- the operation is only carried out by professionally trained personnel.
- the device is used in accordance with these operating instructions.
- extensions, new settings, changes or repairs are only performed by authorised persons.

In the interest of continuous product improvement, we reserve the right to make technical modifications without prior notice.

13. Decommissioning of the product

The DIRECTIVE 2012/19/EU OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 4 July 2012 on waste electrical and electronic equipment (WEEE) lays down measures to protect the environment and human health. Waste electrical and electronic equipment marked as such may not be disposed of as unsorted municipal waste within the European Union. For detailed information regarding the professional disposal of WEEE please contact Dreve or your distributor.

Dreve Dentamid GmbH
Max-Planck-Straße 31 • 59423 Unna/Germany
Tel.: +49 2303 8807-0 • Fax: +49 2303 8807-55
dentamid@dreve.de • dentamid.dreve.de